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ab Seite 3

## **Instruction manual**

starting on page 10

## **Mode d'emploi**

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## **Handleiding**

vanaf pagina 24

DE

GB

FR

NL

## Inhalt

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## Content

Meaning of symbols in these instructions _____	10
Safety instructions _____	11
Proper use _____	11
Risks of injury _____	11
Information for your health _____	11
Notes on use _____	11
Warranty conditions _____	12
Product contents and device overview _____	13
Main unit _____	13
Remote control _____	13
Before first use _____	14
Use _____	14
Cleaning and storage _____	15
Troubleshooting _____	15
Technical data _____	15
Disposal _____	15

---

Dear customer,

We are delighted that you have decided on purchasing our chi massage device. Traditional chi massage activates life force energy (chi). This device is designed to help you harmonise body and spirit to allow your chi to flow freely and if used regularly, you'll feel relaxed yet completely re-energised. We hope you enjoy it!

If you have any questions, please contact customer service via our website:

**[www.service-shopping.de](http://www.service-shopping.de)**

Please carefully read through the operating instructions prior to using the device for the first time and store these instructions in a safe place. These instructions are to accompany the device when it is passed on to others. The manufacturer and importer assume no liability in the event the data in these instructions have not been observed!

As part of ongoing development, we reserve the right to alter products, packaging or enclosed documentation at any time.

## Meaning of symbols in these instructions

All safety notices are marked with this symbol. Please carefully read through and obey the safety notices in order to avoid injury to persons and damage to property.

Tips and recommendations are marked with this symbol.

## Safety instructions

### Proper use

- ☐ The device is not suitable for therapeutic or medical applications.
- ☐ The device is intended for domestic, not commercial use.
- ☐ Only use the device as described in these instructions. Any other use is deemed improper.

### Risks of injury

- ☐ Keep children and animals away from the device and packaging. There is a danger of injury and suffocation!
- ☐ This device is not suitable for use by persons (including children) with restricted sensory or intellectual abilities or with a lack of experience and/or knowledge unless they are supervised by a person responsible for their safety or have received instructions from that person as to how to use the device. Children must be under supervision to ensure that they do not play with the device.
- ☐ Do not use the device if it is damaged.

### Information for your health

- ☐ If you have concerns relating to your health, consult your physician before using the device.
- ☐ Pregnant persons and women should not use the device and/or only after having consulted their physician after having a miscarriage or abortion.
- ☐ Persons with a pacemaker may only use the device after having consulted their physician.
- ☐ Persons suffering from conditions of the heart, kidneys, infarct, high/low blood pressure, anaemia etc. should consult their physician before using the device. If you experience pain when using the device, please stop the exercises and consult your physician.
- ☐ The device is not suitable for persons, who suffer from epilepsy, have unhealed fractures or a disturbance of equilibrium.
- ☐ Patients, who suffer from anaemia or low blood pressure, should place a thin pillow under their head to avoid dizziness.
- ☐ Prior to every work-out, perform some exercises without the device in order to loosen up. Start the exercises slowly and then increase their intensity with time. We expressly point out that using the device improperly or exercising too intensively could be hazardous to your health.

## Notes on use

- ☐ Check the device for damages prior to each use. If the device, the cord or the plug exhibit visible damages, the device may not be used. Do not use the device if it malfunctions, has fallen down or in water.
- ☐ If the device is defective, do not attempt to repair it or replace parts yourself. In the event of damage, contact a professional or customer service.
- ☐ Only connect the device to a properly installed socket with protective earth contacts. The socket must also be easily accessible after the device has been connected. The mains voltage must correspond to the technical data of the device. Only use proper extension cords, whose technical data correspond to those of the device.
- ☐ Only use the device in closed spaces.
- ☐ Do not sit or stand on the vibrating foot rest.
- ☐ Do not insert any objects in the gap between the foot rest and the main unit.
- ☐ Keep body parts, hair etc. away from the openings between the foot rest and main unit.
- ☐ Only use and store the device out of the reach of children and animals. The mains cord may also not be accessible for children. They could pull the device down with it. Do not use the device when water or other liquids are nearby in order to avoid the risk of electrical shock.
- ☐ Maintain adequate distance to heat sources such as ovens or fireplaces to prevent damage to the device.
- ☐ Ensure that the connected mains cord does not pose a tripping hazard.
- ☐ Lay the cable such that it is not crushed or creased and such that it does not come into contact with hot surfaces, including those of the device.
- ☐ In the event of damages/malfunctions, immediately deactivate the device.
- ☐ Do not leave the device on when it is unattended.
- ☐ Only use original spare parts provided by the manufacturer.
- ☐ Turn the device off and pull the mains plug out of the socket
  - if you are not using the device,
  - if a fault occurs during operation,
  - before a storm,
  - before you clean it.
- ☐ If you wish to remove the mains plug from the socket, never pull on the cable, but always on the plug.
- ☐ Do not pull or carry the device by the cord.
- ☐ Never submerge the device in water or other liquids! Ensure that the device, cable or plug do not fall in water or are at risk of getting wet. If the device falls in water, immediately disconnect the power supply. Do not attempt to pull the device out of the water when it is connected to the power supply!
- ☐ Never touch the device, the cable or the plug with wet hands.
- ☐ Do not allow the device to fall and do not subject it to strong percussions.
- ☐ Do not subject the device to extreme temperatures, severe temperature fluctuations or direct sunlight.

## Warranty conditions

Defects arising due to improper handling, damage or attempts at repair are excluded from the warranty. This also applies to normal wear and tear.


## Product contents and device overview

### Main unit

- 1 Vibrating foot rest
- 2 Main unit
- 3 Remote control

### Remote control

- 1 Time control LEDs: Indicate the set time (5, 10, 15 or 20 minutes)
- 2 Speed control LEDs: Indicate the set speed setting
- 3 **Zeit (Time)** button: Button for setting the time interval
- 4 **Geschwindigkeit (Speed)** button: Button for selecting the different speeds
- 5 **An/Aus (On / Off)** button: Button for activating and deactivating the device

## Before first use

### ATTENTION!

- ☐ Keep children and animals away from the packaging. Danger of suffocation exists!

Unpack all parts and make sure that the product contents are complete (see chapter “Product contents and device overview”) and check for any damages from transport. If you identify damages to the components, do not use them (!), instead contact our customer service centre.

## Use

### ATTENTION!

- ☐ If you have concerns relating to your health, consult your physician before using the device.
  - ☐ Discontinue use if you experience pain or discomfort.
  - ☐ Ensure that your clothing doesn't catch in the opening between the foot rest and main unit when in use. Wear comfortable clothing this is not too loose-fitting when using the device.
  - ☐ Keep body parts and hair away from the openings between the foot rest and main unit.
  - ☐ There should be at least 60 minutes between use and your last meal.
- 
1. Set up the device on a stable, level surface. Protect sensitive floors by placing a suitable base under the device.
  2. Only connect the device to a properly installed socket, which corresponds to the device's technical data (also see chapter “Technical data”). Ensure that you lay the mains cable so that it does not pose a tripping hazard.
  3. Lie flat on your back on the floor and place your feet or ankles in the recesses on the vibrating foot rest.
  4. Push the **An/Aus (On/ Off)** button to activate the device. This will immediately activate the device.
  5. Select a speed by repeatedly pushing the **Geschwindigkeit (Speed)** button until you have found a speed that is pleasant for you. The speed control lights indicate which speed you have selected.
  6. Set the timer using the **Zeit (Time)** button. The time control LEDs indicate the time you have set (5, 10, 15 or 20 minutes). Once the preset time has expired, the device automatically deactivates.
  7. Press the **An/Aus (On/ Off)** button to abort use before the set time has elapsed.
  8. After use, pull the mains plug out of the socket.

## Cleaning and storage

### ATTENTION!

- ☐ Prior to each cleaning, pull the plug out of the socket.
  - ☐ Attention! Never submerge the device, the mains cable, and the plug in water or other liquids.
  - ☐ Do not spray any cleaners directly on the device.
  - ☐ When cleaning, do not use any scouring or abrasive cleaners. They could damage the surface.
- 
- If applicable, clean the device with a damp cloth and then dry it off thoroughly.
  - Store the device in a dry area that is out of the reach of children and animals and that is not exposed to direct sunlight or moisture.

## Troubleshooting

If the device is not working properly, please check whether you are able to rectify the problem yourself first. Do not attempt to repair a defective device yourself!

Problem	Possible cause/solution
The device is not working.	<ul style="list-style-type: none"><li>• Is the device activated? Push the On/Off button.</li><li>• Is the mains plug properly connected with the socket?</li><li>• Is the socket defective? Try another socket.</li><li>• Check the fuse for your mains connection.</li></ul>

## Technical data

Model number:	HM01-08Y
Product number:	Z 02937
Voltage supply:	220 – 240 V~ 50/60 Hz
Output:	50 W
Safety class:	I

## Disposal

The packaging material is recyclable. Dispose of the packaging in an environmentally-friendly manner and make it available to the collection service for recyclable materials.

Dispose of the device in an environmentally-friendly manner. It does not belong in the household garbage. Dispose of it at a recycling centre for old electrical and electronic equipment. For more information, please contact the administration in your community.

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All rights reserved.
DS Produkte GmbH, Am Heisterbusch 1, 19258 Gallin, Duitsland